

Rev0-2020.1-Petalinux Project Setup Guide

Author(s): R. Payne

Release: 1.0

Release Date: 21/08/2020

Document Type: Guide

Project: NuPRISM

Keywords: Vivado, Petalinux, NuPRISM, Firmware, Guide, 2020.1

History of Changes

Release No.	Date	Description	Authour(s)
1.0	21/08/2020	Initial Release	R. Payne

Introduction

Vivado and Vivado Project Setup

- 1) If not installed, download Vivado 2020.1 from <https://www.xilinx.com/support/download.html>
 - a. The only chip required is the ZynqMPSoC. To limit the download size, it is recommended that only that one is downloaded
- 2) Clone the project located here: <https://edev-group.triumf.ca/fw/exp/nuprism/xu1/rev0-2020.1-petalinux>
- 3) Open the project using Vivado 2020.1
- 4) Run synthesis
- 5) If you want to set up the Integrated Logic Analyzer (ILA)
 - a. Open synthesized design
 - b. Navigate to set up a debug and select the signals you would like to view
- 6) Run implementation
- 7) Generate Bitstream
- 8) Go to File -> Export -> Export Hardware
 - a. Select Fixed platform type
 - b. Select Include Bitstream
 - c. Export to default border
 - d. Finish

Hardware Setup

Required Hardware:

- Mainboard
- Enclustra XU1
- Mainboard Power Cable
- Xilinx JTAG Debugger
- Power Supply
- Signal Generator
- BNC to Mainboard Cable
- USB-to-Serial Converter
- Ethernet Cable

- 1) Ensure the XU1 is properly in place on the mainboard
- 2) Plug the JTAG Debugger into the Xilinx JTAG port of the mainboard
- 3) Connect your USB-to-Serial Converter to the TX, RX, and GND of the J13 connector
- 4) Connect the two BNC-to-Mainboard Cables to and of J1-4, J6-8, J14-20
- 5) Use your DC power supply to deliver 12V to the mainboard

Insert pic of setup

Petalinux Setup

- 1) Install Ubuntu 18.04.04.
 - a. If using a VM allocate 100GB of memory and 8GB of RAM.
- 2) To download all necessary packages, run:

```
sudo apt-get install -y iproute2 gcc g++ net-tools libncurses5-dev zlib1g:i386 libssl-dev flex bison  
libselinux1 xterm autoconf libtool texinfo zlib1g-dev gcc-multilib build-essential screen pax gawk  
python3 python3-pexpect python3-pip python3-git python3-jinja2 xz-utils debianutils iputils-ping  
libegl1-mesa libsdl1.2-dev pylint3 cpio vim
```
- 3) PetaLinux tools require that your host system /bin/sh is 'bash'. If you are using Ubuntu distribution and your /bin/sh is 'dash', consult your system administrator to change your default system shell /bin/sh with the `sudo dpkg-reconfigure dash` command.
- 4) Download the PetaLinux Tools - Installer - 2020.1 located here:
<https://www.xilinx.com/support/download/index.html/content/xilinx/en/downloadNav/embedded-design-tools.html>
- 5) Run: `mkdir -p ~/pkg/petalinux/pl2020-1`
- 6) In the download location, run: `./petalinux-v2020.1-final-installer.run --dir ../pkg/petalinux/pl2020-1`
- 7) Alter `~/.bashrc` to append: `source ~/pkg/petalinux/pl2020-1/settings.sh`
- 8) Close your terminal and re-open it to see if the tools have installed properly.
- 9) Run the following commands:
 - a. `mkdir ~/project`
 - b. `cd ~/project`
 - c. `petalinux-create -t project -n nuprism --template zynqMP`
 - d. `cd ~/project/nuprism`
 - e. `petalinux-config --get-hw-description <path-to-xsa-exported-from-vivado>`
 - i. Exit the menu that opens, to read more about this see the Petalinux 2020.1 guide
 - f. `petalinux-build`
 - i. This should generate files inside `~/project/nuprism/images/linux/`
 - g. `petalinux-package --boot --u-boot --fpga images/linux/system.bit --format BIN`

- 10) Copy BOOT.BIN, image.ub, and boot.scr to the root partition of your SD card
 - a. To format your SD card, complete the tutorial in Appendix H of the Petalinux 2020.1 user guide:
https://www.xilinx.com/support/documentation/sw_manuals/xilinx2020_1/ug1144-petalinux-tools-reference-guide.pdf
- 11) Extract rootfs.tar.gz to the boot partition of you SD card using: tar xvf rootfs.tar.gz
- 12) Place the SD card into the NuPRISM board, open up a serial connection, and power on the board

MIDAS Host Setup

- 1) Alter ~/.bashrc to append:

```
export GIT_EDITOR="emacs -nw"
export MIDASSYS=$HOME/packages/midas
export MIDAS_EXPTAB=$HOME/online/exptab
export MIDAS_EXPT_NAME=e777
export PATH=$PATH:$MIDASSYS/bin
```
- 2) Run
 - a. mkdir \$HOME/packages
 - b. cd \$HOME/packages
 - c. git clone https://bitbucket.org/tmidas/midas --recursive
 - d. cd midas
 - e. mkdir build
 - f. cd build
 - g. sudo cmake ../
 - h. make install
 - i. cd \$HOME/packages
 - j. git clone https://bitbucket.org/tmidas/rootana
 - k. cd rootana
 - l. git checkout tags/rootana-2020-03-a
 - m. make
 - n. cd \$HOME/packages
 - o. git clone <https://github.com/linev/jsroot.git>
 - p. mkdir \$HOME/online
 - q. cd \$HOME/online
- 3) Alter ~/.bashrc to append:

ROOTANASYS=\$HOME/packages/rootana

JSROOTSYS=\$HOME/packages/jsroot

- 4) Create a file in ~/online/ called exptab with the following contents
 - a. e777 /home/<account>/online <account>
- 5) Run the following 3 programs:
 - a. odbedit
 - b. mhttpd
 - c. mlogger
- 6) You should now be able to open up a browser and navigate to <http://localhost:8080> (or whichever port number is provided by mhttpd) and access the MIDAS frontend panel.

Adding MIDAS to Petalinux and Running the MIDAS Demo

- 1) In the Petalinux project run:
 - a. petalinux-create -t apps --template c++ --name midas --enable
 - b. petalinux-create -t apps --template c++ --name mfe --enable
- 2) Copy the contents of the ds_midas folder into the relative projects generated in project-spec/meta-user/recipe-apps/
- 3) Run:
 - a. petalinux-build -c midas
 - b. petalinux-build -c mfe
- 4) Build the Petalinux project and copy the files over to your SD card as described in the Petalinux setup section (steps 10f,g – 12)
- 5) On you host run
 - a. odbedit
 - b. mhttpd
 - c. mlogger
 - d. mserver
- 6) Once the firmware has booted, run (on the mainboard):
 - a. mfe -h <host-ip-address>:<mserver-address> -e e777
 - i. example: mfe -h 192.168.1.80:1175 -e e777
 - ii. you may need to add the mainboards ip address to a list of trusted ip addresses in the MIDAS control panel
- 7) Start a run by pressing the Start button on the frontend panel start page
- 8) Stop the run and check to see that data has been acquired